




## ASPIRE RENOIR QUATTRO LIFT RECLINER CHAIR

**\$2,671.00**



**SKU:** CHP2275-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Aspire Renoir Quattro Lift Recliner Chair encapsulates modern design, unparalleled comfort, and advanced technology. Created to cater to a broad spectrum of needs, this chair is the epitome of relaxation. It doesn't just recline; it offers a tailored seating experience with its plethora of adjustable features.

### FEATURES

- **Personalized Comfort:** With independent adjustments for the headrest, back, leg, and lumbar support functions, users can tailor their seating experience to their exact preferences.
- **Ergonomic Excellence:** The chair's backrest is designed ergonomically, reducing pressure on the lower sacral area, and ensuring a pain-free, relaxing experience.
- **Spacious Storage:** The large side pockets are perfect for storing reading materials like books and newspapers, keeping them within arm's reach.
- **Universal Design:** The integrated footrest extension accommodates users of various heights, ensuring comfort for everyone.
- **Easy Mobility:** The lift and recline functions simplify the sit-to-stand process, while the 4 swivel castors allow for easy movement of the chair.
- **Tech-Savvy Features:** The backlit hand control with a built-in USB charging point means users can charge their devices without moving. Plus, the chair ensures functionality during power outages with its Lithium Ion Battery Backup.
- **Enhanced Experience:** The adjustable headrest is perfect for those times when you want to watch TV while reclined, ensuring a clear view without compromising comfort.
- **Protective Features:** Both headrest and armrest protectors are included, ensuring longevity.

and cleanliness. Additionally, the handy handset holder ensures the controls are always within reach.

The Aspire Renoir Quattro Lift Recliner Chair is more than just a piece of furniture; it's an experience. With its blend of innovative technology, ergonomic design, and user-focused features, it offers unparalleled comfort and convenience. Whether you're reading, watching TV, or just relaxing, this chair is sure to elevate your relaxation experience. Choose the Renoir Quattro for a harmonious blend of comfort, innovation, and style.

## SPECIFICATIONS

Code	Display Name	Total Width	Total Height	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP227575	Aspire - Renoir Quattro Lift Recliner Chair - Small	750 mm	1010 mm	130 Kg	Dusk Microsuede	510 mm	510 mm	450 mm	450 mm	660 mm
CHP227580	Aspire - Renoir Quattro Lift Recliner Chair - Medium	800 mm	1060 mm	130 Kg	Dusk Microsuede	530 mm	550 mm	470 mm	470 mm	700 mm
CHP227585	Aspire - Renoir Quattro Lift Recliner Chair - Large	830 mm	1060 mm	130 Kg	Dusk Microsuede	530 mm	590 mm	470 mm	470 mm	700 mm