




ASPIRE SHELL CHAIR



SKU: CHP19891-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Shell Chair is a comfortable mobile seating solution designed to maximise user access and social interaction within an Aged Care Facility or Healthcare environment. The compact design allows access under tables for sharing meals, games and other activities. The lightweight frame with anterior tilt function for sit-to-stand transfers reduces the manual handling burden for caregivers.

FEATURES

- Forward tilt functionality assists with chair egress and standing transfers
- Backwards tilt functionality assists in pressure redistribution and patient stability during transfers
- Comfort & Support - Ergonomic wing headrest and lumbar cushioning for improved seating posture
- Storage Pocket - Conveniently carries everyday items such as clipboards, magazines or diaries
- Lightweight and Manoeuvrable - 12.5" rear wheels with attendant foot brake, 5" front swivel castors together with a lightweight design, makes for easy and safe handling for caregivers and a smooth and comfortable ride for users

SPECIFICATIONS

Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
Aspire Shell Chair - 450mm Wide	440 mm	450 mm	450 - 530 mm	135 kg	30 kg

Aspire Shell Chair - 510mm Wide	470 mm	510 mm	450 - 530 mm	135 kg	30 kg
------------------------------------	--------	--------	--------------	--------	-------

BROCHURE

[Aspire Shell Chair Brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au