



## **ASPIRE SHELL CHAIR**

\$3,115.00 Original price was:

\$3,115.00.<u>\$2,804.00</u>Current price is: \$2,804.00.



SKU: CHP19891-H

## **PRODUCT DESCRIPTION**

## **PURPOSE & BENEFITS**

The Aspire Shell Chair is a comfortable mobile seating solution designed to maximise user access and social interaction within an Aged Care Facility or Healthcare environment. The compact design allows access under tables for sharing meals, games and other activities. The lightweight frame with anterior tilt function for sit-to-stand transfers reduces the manual handling burden for caregivers.

#### **FEATURES**

- Forward tilt functionality assists with chair egress and standing transfers
- Backwards tilt functionality assists in pressure redistribution and patient stability during transfers
- Comfort & Support Ergonomic wing headrest and lumbar cushioning for improved seating posture
- Storage Pocket Conveniently carries everyday items such as clipboards, magazines or diaries
- Lightweight and Manoeuvrable 12.5" rear wheels with attendant foot brake, 5" front swivel castors together with a lightweight design, makes for easy and safe handling for caregivers and a smooth and comfortable ride for users

### **SPECIFICATIONS**

Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
Aspire Shell Chair - 450mm Wide	440 mm	450 mm	450 - 530 mm	135 kg	30 kg



Aspire Shell Chair - 510mm Wide	470 mm	510 mm	450 - 530 mm	135 kg	30 kg	
------------------------------------	--------	--------	--------------	--------	-------	--

# **BROCHURE**

**Aspire Shell Chair Brochure**