



ASPIRE SHOWER CHAIR

\$132.00 - \$205.00 Price range: \$132.00 through \$205.00



SKU: BTS1180-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Experience the utmost in safety, comfort, and stability with the Aspire Shower Chair. Our innovative design was crafted to keep user needs at the forefront, resulting in a chair that offers security, adaptability, and a premium shower experience.

FEATURES

- **Ergonomic Contoured Design:** Made from top-grade plastic, the Aspire Chair boasts a contoured seat for maximum postural support, ensuring users remain comfortable and secure throughout their showers.
- **Drainage Holes and Windows:** Strategically positioned for more than just visual appeal. These features ensure rapid water drainage and significantly reduce skin adhesion, offering a dryer and more pleasant seating experience.
- **Height Adjustable Legs:** Personalize your shower experience. A simple button push adjusts the seat height, accommodating varying user needs with ease.
- **Unique Incremental Leg Design:** Uneven shower floors are no match for the Aspire Chair! Thanks to the additional height increments on one leg, the chair remains stable even on uneven surfaces.
- **Corrosion Resistant Frame:** Durability is a promise. Each frame is treated to resist corrosion, ensuring longevity even in the persistent dampness of a shower environment.
- **Raised Armrests:** Transitions made effortless. The elevated armrests are crafted to aid users in transferring on and off the chair seamlessly.
- **Durable Non-slip Feet:** Wet conditions demand enhanced safety. The Aspire Chair feet not only guarantee an excellent grip but also provide unparalleled stability. With an embedded

stainless steel washer, durability is further ensured.

Elevate your shower routine with the Aspire Shower Chair - an epitome of functionality seamlessly blending with comfort. Ideal for both personal and special needs.

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BTS118000	Shower Chair with Arms - Treated Steel	400 mm	440 mm	375 - 550 mm	200 kg	3.4 kg
BTS118010	Shower Chair with Arms - Wide - Treated Steel	400 mm	440 mm	375 - 550 mm	200 kg	3.9 kg