



ASPIRE SHOWER CHAIR



SKU: BTS1180-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Experience the utmost in safety, comfort, and stability with the Aspire Shower Chair. Our innovative design was crafted to keep user needs at the forefront, resulting in a chair that offers security, adaptability, and a premium shower experience.

FEATURES

- **Ergonomic Contoured Design:** Made from top-grade plastic, the Aspire Chair boasts a contoured seat for maximum postural support, ensuring users remain comfortable and secure throughout their showers.
- **Drainage Holes and Windows:** Strategically positioned for more than just visual appeal. These features ensure rapid water drainage and significantly reduce skin adhesion, offering a dryer and more pleasant seating experience.
- **Height Adjustable Legs:** Personalize your shower experience. A simple button push adjusts the seat height, accommodating varying user needs with ease.
- **Unique Incremental Leg Design:** Uneven shower floors are no match for the Aspire Chair! Thanks to the additional height increments on one leg, the chair remains stable even on uneven surfaces.
- **Corrosion Resistant Frame:** Durability is a promise. Each frame is treated to resist corrosion, ensuring longevity even in the persistent dampness of a shower environment.
- **Raised Armrests:** Transitions made effortless. The elevated armrests are crafted to aid users in transferring on and off the chair seamlessly.
- **Durable Non-slip Feet:** Wet conditions demand enhanced safety. The Aspire Chair feet not only guarantee an excellent grip but also provide unparalleled stability. With an embedded

stainless steel washer, durability is further ensured.

Elevate your shower routine with the Aspire Shower Chair - an epitome of functionality seamlessly blending with comfort. Ideal for both personal and special needs.

SPECIFICATIONS

Material	Code	Overall Height (mm)	Overall Width (mm)	Overall Depth (mm)	Seat Height (mm)	Seat Width (mm)	Seat Depth (mm)	Between Arms (mm)	Product Weight (kg)
Steel	BTS118000	735 - 910	512	570	375 - 550	440	400	470	3.4
Wide - Steel	BTS118010	735 - 910	570	570	375 - 550	440	400	525	3.9

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au

