




ASPIRE STAR MOVER - PATIENT MOVER



SKU: LSS390580

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Star Mover offers a safe and comfortable solution for patients and carers doing sit-to-stand transfers. The compact and versatile design is perfect for Homecare, Institutional, and Aged Care environments. A wide seat entry suits a larger range of body shapes. Low footplate height makes patient entry and exit easier.

Large 100 mm rear locking castors have low rolling resistance and offer a smooth ride during transportation.

A narrow base improves access in tight spaces. Height adjustable knee support ensures correct positioning for a larger patient height range. Easy access seat pad handles provide reach-around access to carers. A soft-touch grab bar provides comfortable patient support for improved stability and confidence during transfers.

Width adjustable legs with tapered close position allow access between and, or around chair legs.

FEATURES

- The Aspire Star Mover offers a safe and comfortable solution for patients and carers doing sit to stand transfers.
- The compact and versatile design, is perfect for Homecare, Institutional and Aged Care environments.
- Wide seat entry suits a larger range of body shapes.
- Low footplate height makes patient entry and exit easier.
- Large 100 mm rear locking castors have low rolling resistance and offer a smooth ride during

transportation.

- Narrow base improves access in tight spaces.
- Height adjustable knee support ensures correct positioning for a larger patient height range.
- Easy access seat pad handles provide reach-around access to carers.
- Soft touch grab bar provides a comfortable patient support for improved stability and confidence during transfers.

SPECIFICATIONS

- Length: 925mm
- Height: 1030mm
- Width: 615mm - 900mm
- Knee Pad Width: 390mm
- Knee Pad Height: 260mm
- Knee Pad Height from Base Plate: 250-495mm
- Handle Height: 1030mm
- Seat Pad Height: 685mm - 850mm
- Seat width 500mm - 550mm (Swing back)
- Weight Capacity: 185kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au