




ASPIRE TRANSIT 3 TRANSIT WHEELCHAIR



SKU: MWS452320-H

PRODUCT DESCRIPTION

Purpose & Benefits of the Aspire Transit 3 Transit Wheelchair

The Aspire Transit 3 Transit Wheelchair is a lightweight, durable, and compact mobility solution designed for safe and efficient transportation of individuals with limited mobility. Its folding frame and practical features make it an excellent choice for caregivers and users seeking convenience, reliability, and comfort. Ideal for indoor and outdoor use, this wheelchair provides exceptional manoeuvrability and ease of handling.

Features

- **Lightweight Aluminium Frame:** Designed for easy handling and transportation, the lightweight frame ensures durability without compromising strength.
- **Folding Design:** Compact folding mechanism allows for effortless storage and transport, making it ideal for travel or use in confined spaces.
- **Swing-Away Footrests:** The footrests can swing away or be removed, facilitating easier transfers in and out of the wheelchair.
- **Puncture-Proof Tyres:** Equipped with solid tyres for worry-free use, ensuring smooth operation across various terrains without the risk of punctures.
- **Compact Transit Wheels:** Smaller rear wheels improve manoeuvrability, making it easier to navigate tight indoor spaces.
- **Ergonomic Handgrips with Attendant Brakes:** Features comfortable handgrips and easy-to-use attendant brakes for enhanced safety and control.
- **Padded Upholstery:** The seat and backrest are designed with comfortable, wipe-clean padding for a supportive seating experience.

Specifications

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Frame Weight	Colour	Material
MWS452300	Aspire Transit3 Folding Wheelchair - 400mm	400 mm	40.6cm (16")	440 mm	460 - 500 mm	140 kg	16.4 kg	Silver	Steel
MWS452310	Aspire Transit3 Folding Wheelchair - 450mm	400 mm	45cm (17.7")	440 mm	460 - 500 mm	140 kg	16.6 kg	Silver	Steel
MWS452320	Aspire Transit3 Folding Wheelchair - 500mm	400 mm	50cm (19.6")	440 mm	460 - 500 mm	140 kg	16.8 kg	Silver	Steel

The Aspire Transit 3 Transit Wheelchair is an excellent choice for individuals and caregivers seeking a lightweight, reliable, and easy-to-use wheelchair for everyday mobility needs.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au