




ASPIRE TRI WHEEL WALKER



SKU: WAF705500

PRODUCT DESCRIPTION

Purpose & Benefits of Aspire Tri Wheel Walker

The Aspire Tri Wheel Walker is designed to provide enhanced mobility support for individuals who require assistance with walking, particularly in outdoor environments. This three-wheeled walker offers greater manoeuvrability and stability, making it ideal for users who need a lightweight and versatile mobility aid for navigating both indoor and outdoor spaces.

Features

- **Three-Wheel Design:** The tri-wheel configuration offers superior manoeuvrability and allows for easier navigation around corners and through tight spaces compared to four-wheeled models.
- **Lightweight Frame:** Constructed from lightweight materials, the walker is easy to lift, transport, and handle, reducing the physical effort required by the user.
- **Foldable for Storage:** Features an easy-to-use folding mechanism that allows the walker to be quickly and compactly folded for storage or transport, making it convenient for travel and small living spaces.
- **Adjustable Handles:** Equipped with height-adjustable handles to ensure a comfortable fit for users of various heights, promoting proper posture and reducing strain.
- **Ergonomic Hand Grips:** Includes comfortable, ergonomically designed hand grips that provide a secure hold and reduce hand fatigue during extended use.
- **Storage Bag:** Comes with a detachable storage bag or basket, allowing users to carry personal items, groceries, or other essentials easily, enhancing independence and convenience.
- **Lockable Brakes:** Features lockable hand brakes that provide additional safety and control, allowing the user to securely stop and park the walker when needed.

- **Durable Construction:** Made from high-quality materials that ensure durability and long-lasting performance, suitable for everyday use.

Specifications

- **Safe Working Load:** 115 kg
- **Dimensions:**
 - Height: 825-925 mm
- **Colour Option:** Red
- **Weight:** 5.5 kg

The Aspire Tri Wheel Walker is an excellent choice for individuals seeking a versatile, lightweight, and manoeuvrable mobility aid. Its three-wheel design, combined with practical features such as adjustable handles, lockable brakes, and a convenient storage bag, makes it a valuable tool for enhancing mobility and independence both indoors and outdoors.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au

