




ASPIRE VOGUE ADVENTURE SEAT WALKER

\$807.00



SKU: WAF705450

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Vogue Adventure Seat Walker is a heavy-duty framed walker for users up to 200kg and has been built with Large Multi-Terrain Wheels for durability and stability on and off-road – Front Wheels 11", Back Wheels 9". For increased user safety and stability, Ergonomic locking handbrakes have been added when walking and sitting and comes with a removable utility bag – simple and stylish to remove for shopping – Can hold up to 5kg of essentials. Featuring a wider seat for enhanced support and Compact Turning Circle facilitates improved mobility in tight spaces.

FEATURES

- **Heavy-duty** frame and wider seat provide enhanced support for users up to 200kg.
- **European Design** offers personalised height adjustment, promoting correct posture.
- **Large Multi-Terrain Wheels** for durability and stability on and off road. (Front Wheels 11", Back Wheels 9")
- **Ergonomic Handbrakes** lock for improved user safety and stability.
- **Compact Turning Circle** facilitates improved mobility in tight spaces.
- **Simple Folding Mechanism** enables easy transportation and compact storage (245 mm folded width, perfect for the car).
- **Detachable Storage Bag** carries up to 5kg of essentials.
- **Bonus** walking cane holder provides increased flexibility for your mobility needs.

SPECIFICATIONS

Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight
Aspire Vogue Adventure Walker	770 mm	695 mm	900 - 1035 mm	200 kg	10 kg	Black	230 mm	540 mm	530 mm	530 mm	8.8 kg