



## ASPIRE VOGUE ADVENTURE SEAT WALKER



**SKU:** WAF705450

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Aspire Vogue Adventure Seat Walker is a heavy-duty framed walker for users up to 200kg and has been built with Large Multi-Terrain Wheels for durability and stability on and off-road – Front Wheels 11", Back Wheels 9". For increased user safety and stability, Ergonomic locking handbrakes have been added when walking and sitting and comes with a removable utility bag – simple and stylish to remove for shopping – Can hold up to 5kg of essentials. Featuring a wider seat for enhanced support and Compact Turning Circle facilitates improved mobility in tight spaces.

### FEATURES

- **Heavy-duty** frame and wider seat provide enhanced support for users up to 200kg.
- **European Design** offers personalised height adjustment, promoting correct posture.
- **Large Multi-Terrain Wheels** for durability and stability on and off-road. (Front Wheels 11", Back Wheels 9")
- **Ergonomic Handbrakes** lock for improved user safety and stability.
- **Compact Turning Circle** facilitates improved mobility in tight spaces.
- **Simple Folding Mechanism** enables easy transportation and compact storage (245 mm folded width, perfect for the car).
- **Detachable Storage Bag** carries up to 5kg of essentials.
- **Bonus** walking cane holder provides increased flexibility for your mobility needs.

### SPECIFICATIONS



Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight
Aspire			900 -								
Vogue	770	695	1035	200 kg	10 kg	Black	230 mm	540 mm	530 mm	530 mm	8.8 kg
Adventure Walker	mm	mm	mm								

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)