




ASPIRE VOGUE CARBON FIBRE WALKER

\$873.00



SKU: WAF7053

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Aspire Vogue Carbon Fibre Walker is a sleek, stylish and super lightweight rollator. Its Danish design, profiled tubing and moulded fittings give it a luxurious look and feel. The unique feature of this super lightweight Carbon Fibre walker is that it weighs only 6.2 kg overall and only 4.9 kg frame weight and handles a max user weight of 150 kg.

A key feature of the Vogue walker is its locking mechanism once folded, making it safe and easy to pick up, transport or store. The Aspire Carbon Fibre Seat Walker has ergonomic handles which make holding and controlling the walker comfortable, especially for longer lengths of time.

Adjusting the height of the handles quickly is not a problem with its Auto Locate smooth height adjustment. This feature automatically locates to pre-set heights, making adjusting the height simple and accurate. Being a seat walker, there is also a wide curved backrest for superior comfort and support when you need to sit down. It features a large storage bag that can hold up to 5 kg.

FEATURES

- Smart Posture System equipped with forearm supports that are height, depth, and angle adjustable, the user can achieve optimal centre of gravity for improved balance and weight distribution.
- Luxury look and feel Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours.
- Lightweight and ergonomic design
- “Auto Locate” Smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate.

- The Aspire Vogue Carbon Fibre Walker locks when folded make transportation and storage easier.
- Advanced locking handbrakes for increased user safety when walking and sitting.
- Unbreakable TPE super soft wheels offer high durability and a smooth ride for increased comfort.

If you are looking for a lightweight, stylish, and easy-to-use rollator, the Aspire Vogue Walker is a great option. It is perfect for people who need a little extra support when walking, and it is certain to make your life easier.

SPECIFICATIONS

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight
WAF705350	Seat Walker -											
	Aspire Vogue - Carbon Fibre	680 mm	600 mm	695 - 920 mm	150 kg	6.2 kg	Carbon Fibre	240 mm	460 mm	550 mm	550 mm	4.9 kg
WAF705340	Aspire Vogue Carbon Fibre											
	Seat Walker - Short - Sapphire	655 mm	600 mm	715 - 830 mm	150 kg	5.9 kg	Carbon Sapphire - Blue		460 mm	500 mm	500 mm	4.7 kg
WAF705370	Aspire Vogue Carbon Fibre											
	Seat Walker - Tall - Emerald	655 mm	600 mm	715 - 830 mm	150 kg	6.35 kg	Emerald Green		460 mm	620 mm	620 mm	5.2 kg