



ASPIRE VOGUE FOREARM SEAT WALKER

\$757.00



SKU: WAF705400

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Vogue Forearm Walker offers all the strength, design and integrity of a regular wheeled walker with the added versatility of a compact folding design that allows it to be easily transported or neatly packed away when not in use.

Designed for comfort and usability: Smart Posture System equipped with forearm supports that are height, depth and angle adjustable, the user can achieve optimal centre of gravity for improved balance and weight distribution.

Lightweight and ergonomic design: With small turning circles meaning it is easily manoeuvrable including tight spaces and turns.

Extremely personable: With AutoLocate meaning a smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate. Including adjustable forearm movements in 3 ways!

Need a lighter weight or taller Aspire Walker? [Click here to view our Aspire Vogue Carbon Fibre Walker \(Reg or Tall\)](#)

FEATURES

- The stylish European designed Aspire Vogue Walkers are a lightweight, versatile range with a sturdy frame and easy to fold locking system, suitable for compact transportation when folded.
- Luxury Look and Feel - European design, profiled tubing, moulded fittings and chic colour palette.

- Smart Posture System - Height, angle and depth adjustable forearm supports achieve optimal centre of gravity for improved balance and weight redistribution.
- Smooth height adjustment with AUTOLOCATE - Graduated and numbered upright shaft automatically locates to preset heights, making adjustment simple and accurate.
- Ergonomic Design - Intelligent location of cross brace and generous under seat clearance provide scope for improving patient stride length.
- Rehab Focus - Easy to push with a small turning circle and low rolling resistance to accelerate standing mobility.
- Locks when Folded - Makes lifting the frame easier and minimises storage space required.
- User Comfort - Super-soft wheels and padded forearm surfaces absorb bumps and jolts creating a smooth, comfortable ride.

SPECIFICATIONS

Safe Working Load	150 kg
Total Height	1020 - 1270 mm
Total Length	660 mm
Total Width	620 mm
Seat To Floor Height	620 mm
Seat Width	460 mm
Weight	9.1 kg