




ASPIRE VOGUE LIGHTWEIGHT 2 SEAT WALKER

\$448.20 - \$571.50 Price range: \$448.20 through \$571.50 Sale Available Online Only!



SKU: WAF705310

PRODUCT DESCRIPTION

This Aspire Vogue Lightweight Walker with seat is light in weight and has an ergonomic design. With small turning circles meaning this Lightweight Aspire Vogue Walker is easily manoeuvrable including tight spaces and turns.

ASPIRE VOGUE WALKER PURPOSE & BENEFITS

Extremely personable: With AutoLocate meaning a smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate. Including adjustable forearm movements in 3 ways! This Aspire walker is one of the best [walkers](#).

KEY FEATURES

- **Luxury look and feel** – Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours.
- **Locks when folded** – makes transportation and storage easier.
- **Advanced locking handbrakes** – for increased user safety when walking and sitting.
- **Unbreakable TPE super soft wheels** – offers high durability and a smooth ride for increased comfort.
- **Removable utility bag** – simple and stylish to remove for shopping. Can hold up to 5 kg.

One of the lightest seat walkers on the market! Our Vogue aluminium Seat Walker epitomises function and style with a super lightweight aluminium frame for a superior walking experience.

European Design Ergonomic design offers personalised height adjustment, promoting correct posture. Compact Turning circle for improved mobility in tight spaces. Simple folding mechanism locks for easy transportation and storage. Handbrakes Lock for user safety and stability.

This modern seat walker has a clean and crisp frame design, making it stand out when going for a walk or browsing the shops!

If you are looking for a walker that is lightweight, this one will be perfect for your needs.

This lightweight walker with seat has all the features of an expensive walker without the hefty price tag. It folds down easily and can be stored in small spaces when not in use. The seat on it makes it easy to take a break from standing or sitting while you get around.

Solution: You won't find another lightweight folding aluminium walker with a seat anywhere else!

Buy this Aspire Vogue Lightweight walker with seat online at Medimart using Afterpay or [NDIS method](#) today!

ASPIRE WALKER SPECIFICATIONS

Safe Working Load	150 kg
Bag Safe Working Load	5 kg
Total Height	785 - 920 mm
Total Length	655 mm
Total Width	620 mm
Seat Depth	240 mm
Seat To Floor Height	540 mm
Seat Width	465 mm
Folded Width	220 mm
Wheel Size	8"
Minimum Functional Weight	5.3 kg
Typical Functional Weight	6.45 kg