

ASPIRE WATERFALL DAY CHAIR

\$748.00



SKU: CHP20807-H



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Experience superior comfort and enhanced support with the Aspire Waterfall Day Chair. Designed keeping in mind the orthopedic needs of its users, this chair goes beyond mere seating. With its high back and tailored features, it caters to those who seek a combination of style, adjustability, and unmatched postural support.

FEATURES

- **Contoured, Waterfall Backrest:** The unique segmented backrest design not only enhances the chair's aesthetics but also provides ergonomic support, catering to individual postural needs. Its adjustability means users can achieve their perfect comfort setting.
- **Generously Padded:** Both the seat and armrests come with ample padding, ensuring that users have a comfortable experience and adequate pressure distribution.
- **Easy Maintenance:** With easily removable seat cushions and backrest overlay, cleaning and maintaining the chair becomes a hassle-free task, ensuring hygiene and durability.
- **Robust yet Lightweight:** The chair's sturdy design allows for a high safe working load of up to 200 kg. Yet, it's lightweight enough to ensure easy maneuverability around the home or facility.
- **Customizable Height:** The chair comes with an easy push-button adjustment feature, enabling users or caregivers to modify the height to suit individual needs and preferences.
- **Stable and Secure:** Equipped with non-slip tips, the chair provides enhanced stability, ensuring user safety during both sitting and transferring.
- **Durable Construction:** The powder-coated steel frame is not only easy to clean but is also robust, ensuring that the chair remains a long-lasting addition to any setting.

The Aspire Waterfall Day Chair isn't just a chair; it's a commitment to comfort, support, and

durability. Every feature, from its contoured backrest to its padded surfaces, has been thoughtfully integrated to provide the user with an unparalleled seating experience. Whether you're recovering from surgery, have orthopedic needs, or simply want a comfortable spot to relax, this chair is the ideal choice.

SPECIFICATIONS

| Code | Display Name | Seat Depth | Seat Width | Backrest Height | Seat To Floor Height | Safe Working Load | Colour |
|-----------|---|------------|------------|-----------------|----------------------|-------------------|----------------|
| CHP208070 | Aspire Waterfall Day Chair - Ink Vinyl | 490 mm | 535 mm | 685 mm | 420 - 580 mm | 200 kg | Ink Vinyl |
| CHP208071 | Aspire Waterfall Day Chair - Latte Vinyl | 490 mm | 535 mm | 685 mm | 420 - 580 mm | 200 kg | Latte Vinyl |
| CHP208075 | Aspire Waterfall Day Chair - Oatmeal Fabric | 490 mm | 535 mm | 685 mm | 420 - 580 mm | 200 kg | Oatmeal Fabric |
| CHP208076 | Aspire Waterfall Day Chair - Mocha Fabric | 490 mm | 535 mm | 685 mm | 420 - 580 mm | 200 kg | Mocha Fabric |