




### ASPIRE XL SEAT WALKER



**SKU:** WAF750020SL

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Aspire XL Seat Walker – Silver represents a breakthrough in the world of mobility aids. With its robust design and ergonomic features, this walker is tailored to deliver unmatched comfort and convenience, especially for larger users. Crafted with consideration for users who experience discomfort, arthritis, or joint pain, the Aspire XL Seat Walker ensures every stroll is comfortable and secure.

### FEATURES

- **Durable Build:** Constructed with a robust aluminium frame, the walker promises longevity and reliability. The easy wipe-clean surfaces make maintenance a breeze, ensuring the walker always looks as good as new.
- **Enhanced Weight Capacity:** Designed especially for larger users, this walker features a wide padded seat and boasts an impressive 180kg Safe Working Load.
- **User-Friendly Brakes:** The ultra-soft hand brakes are equipped with palm balls. These are designed to lock effortlessly, making it an ideal choice for those with arthritis or wrist/hand joint pain.
- **Ergonomic Handles:** The anatomically designed, latex-free handles offer a comfortable grip. They are a boon for users with arthritis, ensuring every stroll is pain-free.
- **Customisable Comfort:** The walker is fitted with an adjustable handle height feature. This ensures users can set it to their preferred height, promoting better posture and comfort.
- **Compact Storage:** The easy folding mechanism means that storing the walker or transporting it is hassle-free. This is particularly useful for those who travel or have limited storage space.
- **Convenient Seating:** A wide padded seat provides a comfortable resting place. The seat can

be flipped up to reveal a discreet pocket for storing small items. An additional storage bag is also supplied for larger essentials.

- **Sturdy Wheels:** Solid, puncture-proof 8" wheels promise smooth navigation. Users can be confident while moving around, knowing the wheels are designed for durability and ease of manoeuvre.

The Aspire XL Seat Walker - Silver stands as an embodiment of comfort, durability, and user-focused design. Whether it's the arthritis-friendly features, the ample storage, or the robust build, every aspect of this walker has been meticulously planned to offer users an unmatched walking experience. It's more than just a mobility aid; it's a promise of freedom, safety, and comfort.

## SPECIFICATIONS

Code	Display Name	Seat Width	Total Height	Seat To Floor Height	Safe Working Load	Unit Weight	Colour
WAF750020SL	Aspire XL Seat Walker - Silver	495 mm	850 - 982 mm	560 mm	180 kg	8 kg	Silver

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)