




ABLE SPACER ANTIBACTERIAL



SKU: 3609120AUS

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The ABLE Spacer™ Antibacterial is embedded with silver ion to protect against 99% of all harmful bacteria (demonstrated by ISO 22196 tests). Common plastic spacers can contain up to 40% bacterial contamination that causes harm to the lungs and even pneumonia.

ANTI-BACTERIAL

ANTI-STATIC

BPA FREE

How it works

When bacteria comes into contact with a silver ion protected surface, the silver ions inhibit the bacteria's ability to multiply, therefore they die.

For use with puffer medication. Single patient use.

FEATURES

- Spacer device removes the need of perfect timing of inhalation
- Assists in drug delivery to the lungs
- Improves inhaler technique and drug delivery
- Suitable for all popular MDI's
- Suitable for adults and children
- Ideal volume of 210ml
- Enables tidal breathing
- New hygiene cap

- Instructions on side of device

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au