



ABLE SPACER ANTIBACTERIAL

\$23.00



SKU: 3609120AUS

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The ABLE Spacer™ Antibacterial is embedded with silver ion to protect against 99% of all harmful bacteria (demonstrated by ISO 22196 tests). Common plastic spacers can contain up to 40% bacterial contamination that causes harm to the lungs and even pneumonia.

ANTI-BACTERIAL

ANTI-STATIC

BPA FREE

How it works

When bacteria comes into contact with a silver ion protected surface, the silver ions inhibit the bacteria's ability to multiply, therefore they die.

For use with puffer medication. Single patient use.

FEATURES

- Spacer device removes the need of perfect timing of inhalation
- Assists in drug delivery to the lungs
- Improves inhaler technique and drug delivery
- Suitable for all popular MDI's
- Suitable for adults and children
- Ideal volume of 210ml
- Enables tidal breathing
- New hygiene cap

- Instructions on side of device