

ASTHMA NEBULISER KIT ADULT



\$24.95 Original price was: \$24.95.<u>\$22.00</u>Current price

is: \$22.00.



SKU: nebkitadult

PRODUCT DESCRIPTION

A nebuliser is a device that assists you to breathe in medication as a mist through the use of a mouthpiece or a mask.

A nebuliser transforms a solution of a drug into a light spray.

KEY FEATURES:

- Includes bowl, mask & tubing
- Fits all nebulisers
- New Vixone nebuliser bowl
- Faster delivery of medication
- Simple to use

The user breathes in the mist from the use of ultrasonic power and oxygen to break up the liquid medicine to produce the dose you need.

Treatment - Asthma

Warnings - Always read the label. Follow the directions for use. If symptoms persist, talk to your healthcare professional.

Product - Type Nebuliser

MEDIMART HELP

We provide all the nebulisers you could need, Buy adult medical asthma products online at MediMart using Afterpay or NDIS method today!