




ASTON SHOWER COMMUNE CHAIR



SKU: 530400200-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Aston Shower Commode by Drive Medical is a height adjustable showering commode with 4 locking casters and a waist container. The Aston can be converted into a Self propelling shower commode with the purchase of a self propelling wheel kit.

FEATURES

- Can be used solely as a commode where access to toilet is difficult
- Also assists in transporting users with restricted mobility from toilet to shower, or vice versa, to limit moving and handling transfers
- Self-propelling rear wheel kit available for easy attachment to enable independent use
- Polyurethane (PU) seat, seat-infill, back and armrests are water-impermeable, easy to clean and comfortable
- Bright blue upholstery provides strong colour contrast for visually impaired users
- Seat is moulded for comfort and warmth and has a front cut-away and rear recess to facilitate personal hygiene
- Lightweight, rustproof aluminium frame for ease of manoeuvring/handling and longevity
- Simple e-clip height adjustment to suit different toilet heights
- Tool-free height and angle adjustable footplates accommodate user support requirements
- Footplates swing in/out to facilitate transfers and are easily removable for ease of storage and transportation
- Armrests flip up for side transfers
- Backrest adjusts tool-free in height to suit user and is easily detached/reassembled for storage and transportation

- Supplied with a 5 litre potty with lid, which is easily removable by sliding from rear of seat
- 4 substantial (12.7cm (5")) braked castors, manoeuvre easily over different floor surfaces and can be locked in any direction for improved safety
- Clears standard and disabled WC's for optimum versatility
- Rust-proof, easy-reach aluminium brakes on self-propelling wheels allow stable positioning and are durable
- Can be cleaned easily in domestic environment or thermally decontaminated in automated washing machine
- Optional Attachment - Drive 24" Self Propelling Wheel Kit - [Available Here](#)

SPECIFICATIONS

	Size	Standard	Heavy Duty
Width		56 cm	66 cm
Height		97 to 102cms	97 to 102cms
Depth		90cm	90cm
Seat Width		46cm	56cm
Seat Height		51 to 56cms	51 to 56cms
Seat Depth		43cm	43cm
Weight Limit (Kg)		160kg	260kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au