



## ATMOS ELEVATING LEGREST FOR WHEELCHAIR

**\$292.00**



**SKU:** 15878015-LHS-H

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Atmos Elevating Legrest for Wheelchair** is an adjustable accessory designed to enhance comfort and support for users of the Atmos Ultra Lite Wheelchair range. Available for either left or right side, it allows controlled leg elevation to improve circulation, reduce pressure, and support postural alignment during use.

### Features

- Angle adjustable design enables tailored leg elevation
- Fitted with height and angle adjustable calf pads for comfort and support
- Available individually for left or right side use
- Lightweight construction designed specifically for the Atmos Ultra Lite Wheelchair range
- DVA approved

### Specifications

- **Compatibility:** Atmos Ultra Lite Wheelchair range

- **Options:** Left or Right legrest (sold separately)