



ATMOS ULTRA-LITE SELF PROPELLED WHEELCHAIR



SKU: 15878010-20-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **ATMOS Ultra-Lite Self-Propel Wheelchair** combines lightweight design and high performance, weighing only 8kg for easy transport and storage. Crash-tested to 136kg, it delivers durability, safety, and comfort with an ergonomic frame, cushioned seating, and user-friendly portability features—ideal for short trips and everyday mobility support.

Features

- Ultra-light frame weighing only 8kg, one of the lightest in its class
- Crash tested to 136kg for maximum safety and reliability
- Reinforced frame design enhances durability and strength
- Compact folding frame and collapsible backrest for easy transport
- Pressure-relieving cushioned seat for posture support and comfort
- Removable and washable seat and backrest cushions
- Dual braking system with push-to-lock and user brakes

- Removable anti-tip wheels and swing-away leg rests for safe transfers
- Ergonomic comfort hand grips with integrated storage pockets

Specifications

- **Seat Width Options:** 16" (40cm), 18" (46cm), 20" (50cm)
- **Weight Capacity:** 136kg
- **Unit Weight:** 8kg (without leg rests and cushion)
- **Total Weight:** 10.5kg (with leg rests and cushion)
- **Dimensions:** Length 101cm | Height 92cm | Seat Depth 40cm | Backrest Height 43cm
- **Folded Dimensions:** 68cm (L) x 28cm (W) x 77cm (H)
- **Tyres:** 6" front PU-filled castors, 12" rear PU-filled Air Ride tyres

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au