



ATMOS ULTRA-LITE TRANSIT WHEELCHAIR

\$1,019.70



SKU: 15878000-20-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **ATMOS Ultra-Lite Transit Wheelchair** delivers exceptional portability, comfort, and safety in a lightweight design weighing under 8kg. Ideal for short trips and patient transport, it combines effortless folding and compact storage with a crash-tested reinforced frame, ergonomic comfort features, and superior usability for both user and caregiver.

Features

- Ultra-light frame weighing under 8kg for easy lifting and transport
- Crash tested to 136kg for maximum safety and durability
- Reinforced frame design enhances strength and stability
- One-touch collapsible frame folds into a compact unit for travel or storage
- Collapsible backrest and folding lock strap keep the chair compact and secure in transit
- Cushioned seat contours to the body for posture support and pressure relief
- Removable and washable seat and backrest cushions for easy maintenance



- Ergonomic comfort hand grips with large rear and mesh armrest storage pockets
- Push-to-lock brakes and removable anti-tip wheels ensure stability and safety
- Swing-away leg rests and flip-up footplates with integrated heel loops for easy access
- High-quality solid Air Ride tyres provide smooth mobility
- Available in 16" and 20" seat widths

Specifications

- Seat Width Options: 16" (40cm), 20" (50cm)
- Weight Capacity: 136kg
- **Unit Weight:** 7.5kg (without leg rests and cushion)
- Total Weight: 9kg (with leg rests and cushion)
- Dimensions: Length 87cm | Height 92cm | Seat Depth 40cm | Backrest Height 43cm
- Folded Dimensions: 65cm (L) x 28cm (W) x 73cm (H)
- Tyres: 6" front PU-filled castors, 12" rear PU-filled Air Ride tyres