




## ATMOS ULTRA-LITE TRANSIT WHEELCHAIR



**SKU:** 15878000-20-H

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **ATMOS Ultra-Lite Transit Wheelchair** delivers exceptional portability, comfort, and safety in a lightweight design weighing under 8kg. Ideal for short trips and patient transport, it combines effortless folding and compact storage with a crash-tested reinforced frame, ergonomic comfort features, and superior usability for both user and caregiver.

### Features

- Ultra-light frame weighing under 8kg for easy lifting and transport
- Crash tested to 136kg for maximum safety and durability
- Reinforced frame design enhances strength and stability
- One-touch collapsible frame folds into a compact unit for travel or storage
- Collapsible backrest and folding lock strap keep the chair compact and secure in transit
- Cushioned seat contours to the body for posture support and pressure relief
- Removable and washable seat and backrest cushions for easy maintenance

- Ergonomic comfort hand grips with large rear and mesh armrest storage pockets
- Push-to-lock brakes and removable anti-tip wheels ensure stability and safety
- Swing-away leg rests and flip-up footplates with integrated heel loops for easy access
- High-quality solid Air Ride tyres provide smooth mobility
- Available in 16", 18" 20" seat widths

## Specifications

- **Weight Capacity:** 136kg
- **Unit Weight:** 7.5kg<sup>1</sup>
- **Total Unit Weight:** 9kg<sup>2</sup>
- **Length:** 87cm
- **Height:** 92cm
- **Front Seat Height:** 48cm
- **Seat Width:** 40cm, 46cm, 50cm
- **Seat Depth:** 40cm
- **Backrest Height:** 43cm
- **Folded Length:** 65cm
- **Folded Width:** 28cm
- **Folded Height:** 73cm
- **Height with Backrest Folded:** 73cm

- **Armrest Height to Seat:** 20cm
- **Height with Footplate Removed:** 62cm
- **Leg Rest Adjustment Range:** 10cm
- **Leg Rest Adjustment Range with Optional Extension:** 20cm
- **Front Castor:** 6" PU Filled
- **Rear Wheel:** 12" PU filled

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)

