




INVACARE REA AZALEA MAX SELF PROPELLED WHEELCHAIR

\$7,159.00



SKU: AZALEA-SP-21-MAX

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Engineered to cater specifically to the needs of plus-sized users, the Invacare Rea Azalea Max Wheelchair is a niche product, which complements the existing Invacare Rea Azalea family elegantly. By incorporating features like tilt-in-space with weight-shifting and the DSS (Dual Stability System) mechanism, it leverages proven elements from the broader Invacare Rea Azalea range. Furthermore, the stability provided by this wheelchair is remarkable, thanks to body length compensation in the back and leg rests. Such a feature proves invaluable for users with limited movement or those experiencing difficulty in changing positions.

Moving onto its design, the wheelchair is crafted with a reinforced frame and high-capacity seating, supporting users up to 180 kg effortlessly. Besides ensuring durability and reliability, this robust construction also prioritises user comfort with a specially designed seat and backrest. Also, the seat width and depth are adjustable, thereby accommodating the varying needs of plus-sized users. Therefore, the Invacare Rea Azalea Max Wheelchair doesn't merely make mobility effortless; it also fosters independence and elevates the quality of life for those it serves.

FEATURES

- The Invacare Rea Azalea Max Wheelchair is equipped with a tilt-in-space feature that integrates weight-shifting and a Dual Stability System (DSS). This combination of features allows the user to adjust their seating position effortlessly while maintaining optimal stability and support.
- Ergonomically shaped push handles are closely positioned for ease of handling. This design promotes superior manoeuvrability, offering both the user and caregivers better control in various environments.
- The lower and upper parts of the backrest are reversible, providing an option for additional

space in either the seat or chest area. This adaptability ensures personalised comfort for each individual user.

- The wheelchair features a footplate with a removable midsection. This design facilitates seamless transfers, making the process of entering and exiting the wheelchair simpler and less straining for users and caregivers.
- The Invacare Rea Azalea Max Wheelchair is crash test approved, verifying its high safety standards and robust design, ensuring peace of mind for users and caregivers alike.

SPECIFICATIONS

- Seat width: 55cm (21.6") – 66cm (25.9")
- Seat depth: 53 – 60 cm
- Seat height: 40 – 45 cm
- Total width: seat width plus 25 cm
- Total length: 112 cm
- Legrest angle: adjustable
- Tilt range: 25°
- Transport weight: 33 kg
- Total weight unloaded: 54kg
- Turning radius: 77 cm
- Maximum safe slope: 7°
- Maximum user weight: 180 kg

BROCHURE

[Invacare Azalea Max Brochure](#)

COMPLIANCE

ARTG 99518

AS/NZS 3695

Crash Tested: ISO 7176-19

WARRANTY

Frame: 2 years | Components: 2 years