




INVACARE REA AZALEA SELF PROPELLED WHEELCHAIR

\$4,790.00 - \$6,950.00 Price range: \$4,790.00 through \$6,950.00



SKU: AZALEA-SP-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Azalea Wheelchair - Self Propelled provides unrivalled flexibility, allowing for tailored fitting to each user's unique needs, promoting ease of mobility and independence. As an ideal companion for daily activities, this wheelchair's design combines comfort, adjustability, and user-friendly operation, positioning it as the optimal choice for those in need of reliable mobility solutions.

FEATURES

- The Azalea Wheelchair - Self Propelled features a Flex 3 tension adjustable backrest and backrest plates, providing superior customisation to ensure comfort and proper posture.
- The 25° seat tilt and 30° backrest recline allow users to find the most comfortable sitting position, reducing strain and optimising long-term usage.
- Unique to this model is the biomechanical positioning hinge on the backrest, promoting natural movement and increased mobility for the user.
- The wheelchair is available in 3 seating width options, offering a personalised fit to accommodate various body types.

SPECIFICATIONS

- Seat Width: 40.6cm (16"); 45.72cm (18"); 51cm (20"); 61cm (24");
- Width: seat width + 23cm
- Length: 102cm
- Product Weight: 34kg
- Warranty: 2 years
- Max User Weight: 135kg, 180kg on 24" Max

- Seat Width: 40.6cm (16")
- Seat Depth: 43cm - 50cm
- Seat Height: 45cm
- Tilt Range: 25°

COMPLIANCE

ARTG 99518

AS/NZS 3695

Crash Tested: ISO 7176-19

WARRANTY

Frame: 2 years | Components: 2 years