



# **BARDWELL HOIST MAT**



SKU: 1001028-H

## PRODUCT DESCRIPTION

#### **PURPOSE AND BENEFITS**

The Hoist mat surface is ultra clear and made from twice-extruded heavy duty PVC, designed for mechanical hoists, walkers and heavy load bearing. The Hoist mat is secured to the floor with rounded grippers or an adhesive backing.

The Hoist mat creates a strong flat stable surface on top of a carpeted floor that enables you to place a mechanical hoist on top and to safely lift patients in and out of hospital beds.

#### **FEATURES**

- The Hoist mat is for bed ridden patients who depend on a hoist to get out of bed or a machine to assist them with walking.
- A smooth clear heavy duty PVC surface with rounded grippers or an adhesive backing on the underside, creates a strong and stable flooring base with the load bearing capacity to support a hoist carrying the weight of a patient.
- Ideal for hospitals, aged care, and patients private rooms.

### **SIZING**

- 1.5m x 1.2m CMDIPLOMAT
- 2.4m x 1.2m HOISTMAT

For further information go to **medimart.com.au** or give us a call: **1300 633 063** 



# SA

# **Brighton Showroom**

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

# **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

# **Reynella Showroom**

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

# **VIC**

#### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

# NT

#### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au

# **Monash Showroom**

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au