



BARIATRIC KING COMFORT HEIGHT ADJUSTABLE SUPPORT CHAIR



SKU: AGE300B

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The King Comfort Bariatric Low Back Chair offers superb seating support with its low back rest and adjustable height. Firm padded armrests assist with transferring in and out of the chair. The King Comfort is suitable for larger clients with its option of seat widths.

FEATURES

- King Comfort Bariatric Healthcare Chair
- Custom Designed Heavy-Duty Metal Frame
- Frames are Finished in a Durable Powder coat
- Seat & Back Rests are Solid Ply Timber
- Foams are Dunlop Enduro Fire Retardant
- Padded Armrests
- Excellent Comfort & Support
- Commercial Grade Upholstery
- High-Quality Stitching
- Hand Built in Australia
- Floor to Seat Adjustment

COLOUR & STYLE

Some colours in the image have been discontinued, select a colour from the available list of options.



LOOKING FOR A DIFFERENT VERSION?

- Standard King Comfort Height Adjustable Support Chair
- Broad King Comfort Height Adjustable Support Chair

DIMENSIONS

	<u>Standard</u>	<u>Broad</u>	Bariatric
Width	450mm	550mm	650mm
Depth	450mm	450mm	450mm
Floor to seat Adjustment	460mm-585mm	460mm-585mm	460mm-585mm
Back Rest Height	370mm	370mm	370mm
Arm Rest Height	180mm	180mm	180mm
Weight Capacity	350kg	350kg	350kg
Package Weight	12kg	12kg	14kg

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom
343 Brighton Rd, North
Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771 reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au



NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au