



## **BARIATRIC KING COMFORT HEIGHT ADJUSTABLE SUPPORT CHAIR**

**\$408.00**



**SKU: AGE300B**

### **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The King Comfort Bariatric Low Back Chair offers superb seating support with its low back rest and adjustable height. Firm padded armrests assist with transferring in and out of the chair. The King Comfort is suitable for larger clients with its option of seat widths.

### **FEATURES**

- King Comfort Bariatric Healthcare Chair
- Custom Designed Heavy-Duty Metal Frame
- Frames are Finished in a Durable Powder coat
- Seat & Back Rests are Solid Ply Timber
- Foams are Dunlop Enduro Fire Retardant
- Padded Armrests
- Excellent Comfort & Support
- Commercial Grade Upholstery
- High-Quality Stitching
- Hand Built in Australia
- Floor to Seat Adjustment

### **COLOUR & STYLE**

Some colours in the image have been discontinued, select a colour from the available list of options.

## LOOKING FOR A DIFFERENT VERSION?

- [Standard King Comfort Height Adjustable Support Chair](#)
- [Broad King Comfort Height Adjustable Support Chair](#)

## DIMENSIONS

|                                 | <a href="#">Standard</a> | <a href="#">Broad</a> | <b>Bariatric</b> |
|---------------------------------|--------------------------|-----------------------|------------------|
| <b>Width</b>                    | 450mm                    | 550mm                 | 650mm            |
| <b>Depth</b>                    | 450mm                    | 450mm                 | 450mm            |
| <b>Floor to seat Adjustment</b> | 460mm-585mm              | 460mm-585mm           | 460mm-585mm      |
| <b>Back Rest Height</b>         | 370mm                    | 370mm                 | 370mm            |
| <b>Arm Rest Height</b>          | 180mm                    | 180mm                 | 180mm            |
| <b>Weight Capacity</b>          | 350kg                    | 350kg                 | 350kg            |
| <b>Package Weight</b>           | 12kg                     | 12kg                  | 14kg             |