



# DRIVE DEVILBISS BARIATRIC STEEL TRANSPORT CHAIR - BLACK 22" - 200KG

\$472.00



**SKU:** BTR22BLKAU

### **PRODUCT DESCRIPTION**

## **Purpose & Benefits**

The **Drive DeVilbiss Bariatric Steel Transport Chair - Black 22**" is built for comfort, strength, and stability. Designed with a dual reinforced steel cross brace and heavy-duty upholstery, it supports users up to 200 kg while providing a safe and durable transport solution.

#### **Features**

- Dual reinforced steel cross braces for added support and stability
- Heavy-duty steel frame construction with a 200 kg weight capacity
- Durable nylon reinforced upholstery with back carry pouch for storage
- Removable desk-length armrests for easy side transfers and access to surfaces
- Height adjustable footrests to accommodate user needs
- 12" rear solid tyres for smooth, low-maintenance mobility
- Stepper tube allows attendants to safely navigate kerbs



# **Specifications**

• **SKU:** BTR22BLKAU

• Seat Width: 56 cm (22")

• Seat Depth: 44 cm (17.3")

• Seat Height: 52 cm (20.5")

• Overall Height: 92 cm (36.2")

• Overall Width: 71 cm (27.9")

• Overall Length: 104 cm (40.9") including footrests

• Width (Closed): 29 cm (11.4")

• Product Weight: 22 kg

• Maximum User Weight: 200 kg (31.5 st)