



DRIVE DEVILBISS BARIATRIC STEEL TRANSPORT CHAIR - BLACK 22" - 200KG



SKU: BTR22BLKAU

PRODUCT DESCRIPTION

Purpose & Benefits

The **Drive DeVilbiss Bariatric Steel Transport Chair - Black 22"** is built for comfort, strength, and stability. Designed with a dual reinforced steel cross brace and heavy-duty upholstery, it supports users up to 200 kg while providing a safe and durable transport solution.

Features

- Dual reinforced steel cross braces for added support and stability
- Heavy-duty steel frame construction with a 200 kg weight capacity
- Durable nylon reinforced upholstery with back carry pouch for storage
- Removable desk-length armrests for easy side transfers and access to surfaces
- Height adjustable footrests to accommodate user needs
- 12" rear solid tyres for smooth, low-maintenance mobility
- Stepper tube allows attendants to safely navigate kerbs

Specifications

- **SKU:** BTR22BLKAU
- **Seat Width:** 56 cm (22")
- **Seat Depth:** 44 cm (17.3")
- **Seat Height:** 52 cm (20.5")
- **Overall Height:** 92 cm (36.2")
- **Overall Width:** 71 cm (27.9")
- **Overall Length:** 104 cm (40.9") including footrests
- **Width (Closed):** 29 cm (11.4")
- **Product Weight:** 22 kg
- **Maximum User Weight:** 200 kg (31.5 st)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au