



## BATHMASTER DELTIS BATH LIFT



**SKU:** PAT-091324169

## PRODUCT DESCRIPTION

### Purpose & Benefits of Bathmaster Deltis Bath Lift

The Bathmaster Deltis Bath Lift is designed to provide a safe and comfortable bathing experience for individuals with mobility limitations, allowing them to easily lower themselves into and out of the bathtub with minimal assistance. This innovative bath lift offers several features that promote convenience, independence, and relaxation during bathing.

### Features

- **Powered Lifting Mechanism:** Utilises a rechargeable battery-powered lifting mechanism to smoothly lower and raise the user in and out of the bathtub, eliminating the need for manual lifting or assistance from caregivers.
- **Compact and Lightweight Design:** Features a compact and lightweight design that makes it easy to transport and store when not in use, ideal for users with limited bathroom space or for travel.
- **Comfortable Seat and Backrest:** Equipped with a comfortable seat and supportive backrest that provide stability and comfort during bathing, reducing strain on the user's body.
- **Adjustable Recline:** Offers adjustable recline positions to accommodate individual comfort preferences and provide optimal support while bathing.

- **Waterproof Construction:** Made from waterproof materials that are resistant to moisture and corrosion, ensuring durability and longevity in wet environments.
- **Easy-to-Use Hand Controller:** Includes a simple-to-operate hand controller with large, ergonomic buttons that allow users to control the lifting and reclining functions with ease.
- **Safety Features:** Incorporates safety features such as a low battery indicator, emergency stop button, and stability locking mechanism to ensure safe and secure operation during bath transfers.

## Specifications

- **Safe Working Load:** 140kg
- **Dimensions:**
  - Backrest: 680mm (H) x 350mm (W), 3.7kg
  - Back angle: Tilts to 40 degrees
  - Base: 530mm (D) x 220mm (W)
  - Seat dimensions: 500mm (D) x 375mm (W)
  - Seat height: 65mm to 460mm
  - Seat weight: 6.8kg
- **Battery Life:** Offers long-lasting battery life to provide multiple bath cycles on a single charge, ensuring uninterrupted operation during bathing.
- **Material:** Constructed from high-quality, durable materials that are easy to clean and maintain for optimal hygiene.
- **Safety Certification:** Meets rigorous safety standards and regulations to ensure user safety and peace of mind.

The Bathmaster Deltis Bath Lift offers a safe, convenient, and relaxing bathing solution for individuals with mobility challenges, allowing them to maintain their independence and dignity while enjoying a soothing bath experience. With its user-friendly design and advanced features, this bath lift enhances safety and comfort in the bathroom, promoting overall well-being and quality of life for users and caregivers alike.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)