



BELLAVITA BATH LIFT

\$1,011.00



SKU: S11499

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Enjoy a bath with confidence using the lightest battery powered bath lift available – the Bellavita bath lift. At the touch of a button using the hand control of this ergonomically-designed bath lift, users with restricted mobility are gently lowered down into the bath. The bath lift also lifts the user to the top again for safely getting out.

Innovative and effective: Made in Germany, this ultra-modern bath lift takes the clinical look out and is more in keeping with today's requirements. Developed with the help of Nursing Professionals, it includes several excellent new features.

Lightweight yet strong: Lightest battery-powered reclining bathlift with a unique front aperture for improved personal hygiene.

Record defining: Ergonomically designed for you to enjoy a great soak in the bath with a record seat height of 48cm, very low seat travel of 6cm and a backrest that reclines to a market leading 50 degrees.

KEY FEATURES

- The backrest reclines to a market-leading 50 degrees to accommodate customers requiring a greater hip-to-back sitting position, such as those with stiff hips.
- The cutout in the seat's front facilitates personal hygiene.
- Can be used with or without covers;

washable hygienic covers are simply fitted/removed and offer additional bathing comfort (available in white or blue)

- Extremely low seat travel of 6cm (2.4") from the bottom of the bath enhances water immersion and is perfect for shallow baths.
- Record seat height of 48cm (18.9") allows for safe level transfers in deeper baths • Weighs only 9.3kg (20.5lb) to make transporting and handling easier.
- Quickly and simply separates into two lightweight halves for simple installation/removal, storage, and travel.
- Contoured, slim-line profile facilitates handling and provides a contemporary, aesthetically pleasant appearance that complements most bathrooms.
- Four quick-release suction feet permit combat-free removal from bath.
- Lightweight, completely waterproof, floatable hand control with easy-to-use buttons.
- The handheld remote's lithium-battery technology allows for convenient recharging in numerous locations.
- An indicator light illuminates when charging is necessary, and a safety mechanism stops the bath lift from lowering without sufficient energy to return it to the top of the bath.
- The emergency stop button built into the handle ensures maximum safety.
- User-friendly jack plug connection to charger for simplicity of usage.
- Produced from high-quality, lightweight plastic that is recyclable and simple to clean.
- Weight limit of 140 kg.

SPECIFICATIONS

| | |
|--|---|
| Max user weight (SWL) | 140kg |
| Total weight (excl. hand control) | 9.2kg |
| Overall length | 83cm (32.5 inches) – with backrest in a reclined position |
| Overall width | 70cm (27.5 inches) – seat including side flaps |
| Seat width | 35cm (14 inches) |
| Seat depth | 42cm (16.5 inches) |
| Seat length | 57cm (22.5 inches) |
| Seat height | 6cm – 48cm (2.5 – 19 inches) |
| Backrest height | 64cm (25 inches) |
| Backrest width | 32 – 36cm (13 – 14 inches) |
| Backrest reclines | 50 degrees |
| Baseplate width | 31cm (12 inches) |
| Baseplate length | 56cm (22 inches) |
| Warranty | 2 years warranty on battery/handset. |
| Colour | Blue |

Buy this Bathlift online at Medimart using Afterpay or [NDIS method](#) today!

