




BELLAVITA BATH LIFT



SKU: S11499

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Enjoy a bath with confidence using the lightest battery powered bath lift available – the Bellavita bath lift. At the touch of a button using the hand control of this ergonomically-designed bath lift, users with restricted mobility are gently lowered down into the bath. The bath lift also lifts the user to the top again for safely getting out.

Innovative and effective: Made in Germany, this ultra-modern bath lift takes the clinical look out and is more in keeping with today's requirements. Developed with the help of Nursing Professionals, it includes several excellent new features.

Lightweight yet strong: Lightest battery-powered reclining bathlift with a unique front aperture for improved personal hygiene.

Record defining: Ergonomically designed for you to enjoy a great soak in the bath with a record seat height of 48cm, very low seat travel of 6cm and a backrest that reclines to a market leading 50 degrees.

KEY FEATURES

- The backrest reclines to a market-leading 50 degrees to accommodate customers requiring a greater hip-to-back sitting position, such as those with stiff hips.
- The cutout in the seat's front facilitates personal hygiene.
- Can be used with or without covers;

washable hygienic covers are simply fitted/removed and offer additional bathing comfort (available in white or blue)

- Extremely low seat travel of 6cm (2.4") from the bottom of the bath enhances water immersion and is perfect for shallow baths.
- Record seat height of 48cm (18.9") allows for safe level transfers in deeper baths • Weighs only 9.3kg (20.5lb) to make transporting and handling easier.
- Quickly and simply separates into two lightweight halves for simple installation/removal, storage, and travel.
- Contoured, slim-line profile facilitates handling and provides a contemporary, aesthetically pleasant appearance that complements most bathrooms.
- Four quick-release suction feet permit combat-free removal from bath.
- Lightweight, completely waterproof, floatable hand control with easy-to-use buttons.
- The handheld remote's lithium-battery technology allows for convenient recharging in numerous locations.
- An indicator light illuminates when charging is necessary, and a safety mechanism stops the bath lift from lowering without sufficient energy to return it to the top of the bath.
- The emergency stop button built into the handle ensures maximum safety.
- User-friendly jack plug connection to charger for simplicity of usage.
- Produced from high-quality, lightweight plastic that is recyclable and simple to clean.
- Weight limit of 140 kg.

SPECIFICATIONS

Max user weight (SWL)	140kg
Total weight (excl. hand control)	9.2kg
Overall length	83cm (32.5 inches) – with backrest in a reclined position
Overall width	70cm (27.5 inches) – seat including side flaps
Seat width	35cm (14 inches)
Seat depth	42cm (16.5 inches)
Seat length	57cm (22.5 inches)
Seat height	6cm – 48cm (2.5 – 19 inches)
Backrest height	64cm (25 inches)
Backrest width	32 – 36cm (13 – 14 inches)
Backrest reclines	50 degrees
Baseplate width	31cm (12 inches)
Baseplate length	56cm (22 inches)
Warranty	2 years warranty on battery/handset.
Colour	Blue

Buy this Bathlift online at Medimart using Afterpay or [NDIS method](#) today!

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au