


BELLMAN ALARM CLOCK PRO



SKU: BE1370-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Alarm Clock Pro** from Bellman & Symfon provides reliable wake-up alerts through sound, vibration, and flashing lights. Designed for individuals who are deaf, hard of hearing, or heavy sleepers, it combines a loud alarm, powerful bed shaker, and high-intensity light signals. With features like a night light, smart snooze, and backup power, it ensures dependable performance every day.

FEATURES

- Extra-loud 100dB alarm that increases gradually in volume
- Frequency-sweeping tone ensures all hearing ranges are reached
- High-intensity LED flashing lights provide visual wake-up alerts
- Strong vibrating bed shaker (included with certain models) emits both vibration and sound
- Night light beacon with gentle blue illumination for easy navigation in the dark
- Smart snooze function reduces intervals from 9 to 2 minutes for gradual wake-up
- Telephone line connection alerts you to incoming calls

- External input socket for connection to Bellman Visit accessories such as a mobile phone sensor or contact mat
- Mains powered with 4 × AAA NiMH backup batteries that maintain full alert function during power loss
- Large LCD display with adjustable backlight and simple push-rotate dials

AVAILABLE VARIANTS

Model	Description
Bellman Alarm Clock Pro	Standard model with loud alarm, flashing lights, and night light.
Bellman Alarm Clock Pro + Bed Shaker	Includes vibrating bed shaker for additional tactile alerts.

SPECIFICATIONS

- **Alarm Volume:** Up to 100dB
- **Alert Types:** Sound, light, vibration (model dependent)
- **Power Source:** Mains powered with 4 × AAA NiMH battery backup
- **Connectivity:** Telephone line and external socket for Bellman Visit accessories
- **Display:** Large LCD with adjustable backlight
- **Night Light:** Gentle blue light beacon for dark environments

Note: Compatible with Bellman Visit system accessories including mobile phone sensor and contact mat (sold separately).

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au