

BELLMAN AUDIO EARBUDS



SKU: BE9124-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Audio Earbuds** are high-quality stereo in-ear headphones designed to deliver clear, natural sound for everyday listening. Ideal for use with Bellman & Symfon devices, they provide comfort, clarity, and lightweight convenience. The version with a built-in microphone is specially designed for the Maxi Pro personal amplifier, allowing you to listen to music or take phone calls seamlessly.

FEATURES

- Comfortable, lightweight stereo in-ear design
- Delivers clear, balanced sound across a wide frequency range
- Available with or without an inline microphone
- Inline microphone version includes answer/end call button for phone use
- Compatible with the Bellman Maxi Pro personal amplifier (microphone model only)
- 3.5mm angled stereo plug for secure and easy connection

AVAILABLE VARIANTS

Model	Description
Bellman Audio Earbuds (No Mic)	Standard stereo earbuds without a microphone – suitable for music and audio playback.
Bellman Audio Earbuds with Inline Microphone (Maxi Pro Only)	Includes integrated microphone and call control button – designed for the Bellman Maxi Pro amplifier.

SPECIFICATIONS

- **Impedance:** 16 Ω
- **Sensitivity:** 93 dB/mW at 1 kHz
- **Frequency Range:** 20 – 20,000 Hz
- **Connector:** 3.5 mm angled stereo plug
- **Weight:** 16 g
- **Colour:** Black
- **Cable Length:** 1.25 m

Note: The inline microphone version is compatible only with the **Bellman Maxi Pro** personal amplifier.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au