

BELLMAN AUDIO HEADPHONES WITH INLINE MICROPHONE (MAXI PRO ONLY)



SKU: BE9233



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Audio Headphones with Inline Microphone (BE9229)** are designed specifically for use with the Bellman Maxi Pro personal amplifier. These comfortable on-ear stereo headphones deliver high-quality sound for listening to music, phone calls, or conversations without distractions. The built-in microphone and answer/end call button make communication simple and hands-free, providing a clear and convenient listening experience for everyday use.

FEATURES

- Optimised for use with the Bellman Maxi Pro personal amplifier
- Integrated microphone with answer/end call button for hands-free operation
- Delivers clear, high-fidelity stereo sound for calls and media
- Comfortable on-ear design suitable for extended wear
- Durable, lightweight build with tangle-resistant cable
- Standard 3.5 mm angled stereo connector for secure connection

SPECIFICATIONS

- **Model:** BE9229
- **Impedance:** $32 \Omega \pm 10\%$
- **Sensitivity:** $113 \text{ dB/mW} \pm 3 \text{ dB at } 1 \text{ kHz}$
- **Frequency Range:** 20 – 20,000 Hz
- **Maximum Power Handling:** 100 mW
- **Connector:** 3.5 mm angled stereo teleplug
- **Driver Unit:** Ø 30 mm anisotropic magnet

- **Cable Length:** 1.2 m
- **Weight:** 63 g
- **Colour:** Black

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au