

BELLMAN AUDIO HIGH POWER NECKLOOP



SKU: BE9159

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Audio High Power Neckloop** is designed for hearing aid users who require strong, clear audio transmission from personal amplifiers, smartphones, computers, or other audio devices. Worn comfortably around the neck, it delivers high-output, dual-channel sound directly to hearing aids set to the “T” (telecoil) position. Ideal for individuals seeking amplified clarity during calls, music, or conversations, this lightweight and hygienic neckloop combines reliability, comfort, and safety.

FEATURES

- High-output dual-channel transmission for powerful, clear sound
- Soft fabric tube design ensures comfort and hygiene during long use
- Breakaway safety strap for enhanced protection and easy release
- Tangle-free cable for convenient daily use
- Gold-plated angled 3.5 mm plug resists corrosion and ensures stable connection
- Compatible with personal amplifiers, smartphones, tablets, and computers

SPECIFICATIONS

- **Model:** BE9159 / BE9161
- **Input Impedance:** $2 \times 5 \Omega$
- **Connector:** 3.5 mm gold-plated 90° angled stereo teleplug
- **Output Power:** 1500 mA/m @ 15 cm distance with 2×50 mW input
- **Dimensions (W × H × D):** $3.5 \times 4 \times 1$ cm
- **Weight:** 60 g
- **Neck Strap Length:** 70 cm

- **Cable Length:** 90 cm (BE9159) / 15 cm (BE9161)
- **Neckloop Diameter:** 22 cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au