



BELLMAN AUDIO HIGH POWER NECKLOOP

\$120.00 Original price was: \$120.00.\$108.00 Current price is: \$108.00.



SKU: BE9159

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Audio High Power Neckloop** is designed for hearing aid users who require strong, clear audio transmission from personal amplifiers, smartphones, computers, or other audio devices. Worn comfortably around the neck, it delivers high-output, dual-channel sound directly to hearing aids set to the "T" (telecoil) position. Ideal for individuals seeking amplified clarity during calls, music, or conversations, this lightweight and hygienic neckloop combines reliability, comfort, and safety.

FEATURES

- High-output dual-channel transmission for powerful, clear sound
- Soft fabric tube design ensures comfort and hygiene during long use
- Breakaway safety strap for enhanced protection and easy release
- Tangle-free cable for convenient daily use
- Gold-plated angled 3.5 mm plug resists corrosion and ensures stable connection
- Compatible with personal amplifiers, smartphones, tablets, and computers

SPECIFICATIONS

Model: BE9159 / BE9161
Input Impedance: 2 × 5 Ω

• Connector: 3.5 mm gold-plated 90° angled stereo teleplug

• Output Power: 1500 mA/m @ 15 cm distance with 2 × 50 mW input

• Dimensions (W \times H \times D): 3.5 \times 4 \times 1 cm

• Weight: 60 g

• Neck Strap Length: 70 cm



• Cable Length: 90 cm (BE9159) / 15 cm (BE9161)

• Neckloop Diameter: 22 cm