

## BELLMAN BED SHAKER



**SKU:** BE1270

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The **Bellman Bed Shaker** (BE1270) provides powerful, reliable vibration alerts to wake even the heaviest sleepers. Designed to work with Bellman & Symfon alarm clocks and Visit receivers, it offers a safe, low-voltage solution for individuals who are deaf, hard of hearing, or deep sleepers. Simply place it under your pillow or mattress to be awakened by strong vibrations when your connected device is triggered.

### FEATURES

- Powerful vibration ensures effective wake-up for deep or hard-of-hearing sleepers
- Anti-slip protection keeps the shaker securely in place during operation
- Safe low-voltage design for worry-free use
- 2-metre cable provides flexible placement under the pillow or mattress
- Compatible with all Bellman Visit receivers and alarm clocks
- No internal battery required — powered directly from connected device

### SPECIFICATIONS

- **Model:** BE1270
- **Operating Voltage:** 2.0 – 4.0 V DC (from Bellman Visit receiver)
- **Power Consumption:** 250 – 750 mA
- **Dimensions (W × H × D):** 88 × 88 × 27 mm
- **Weight:** 120 g
- **Cable Length:** 2 m
- **Colour:** White with red rubber anti-skid protector

- **Usage:** Indoor use only

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriiffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)