

BELLMAN BED SHAKER



\$65.00 Original price was: \$65.00.<u>\$58.50</u>Current price

is: \$58.50.



SKU: BE1270

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Bed Shaker** (BE1270) provides powerful, reliable vibration alerts to wake even the heaviest sleepers. Designed to work with Bellman & Symfon alarm clocks and Visit receivers, it offers a safe, low-voltage solution for individuals who are deaf, hard of hearing, or deep sleepers. Simply place it under your pillow or mattress to be awakened by strong vibrations when your connected device is triggered.

FEATURES

- Powerful vibration ensures effective wake-up for deep or hard-of-hearing sleepers
- Anti-slip protection keeps the shaker securely in place during operation
- Safe low-voltage design for worry-free use
- 2-metre cable provides flexible placement under the pillow or mattress
- Compatible with all Bellman Visit receivers and alarm clocks
- No internal battery required powered directly from connected device

SPECIFICATIONS

• **Model:** BE1270

• Operating Voltage: 2.0 - 4.0 V DC (from Bellman Visit receiver)

• Power Consumption: 250 - 750 mA

• Dimensions (W \times H \times D): 88 \times 88 \times 27 mm

• Weight: 120 g

• Cable Length: 2 m

• Colour: White with red rubber anti-skid protector



• **Usage:** Indoor use only