

BELLMAN LEAD FOR MINO CHARGER



SKU: BE9148

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Lead for Mino Charger** is a genuine replacement charging cable designed for the Bellman Audio Mino personal amplifier. It connects the Mino amplifier to a USB power source, ensuring efficient and reliable charging whether at home or on the go. Ideal as a spare or replacement cable to keep your device powered at all times.

FEATURES

- Official Bellman replacement charging lead for the Mino personal amplifier
- Connects easily to any standard USB power source
- Angled connector for a secure and convenient fit
- Durable, tangle-resistant cable design for everyday use
- Compact and lightweight — perfect for travel or home use

SPECIFICATIONS

- **Compatibility:** Bellman Audio Mino personal amplifier
- **Connector Type:** USB-A to DC barrel plug (angled)
- **Cable Length:** Approximately 1 metre
- **Colour:** Black
- **Usage:** For charging only — not for data transfer

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au