


## **BELLMAN VIBIO - WIRELESS BLUETOOTH BED SHAKER ALARM**



**SKU:** BE1221

### **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The **Bellman Vibio Wireless Bluetooth Bed Shaker Alarm** is a smart, portable solution for reliable wake-ups, designed especially for deaf users, those with hearing loss, and heavy sleepers. Connecting via Bluetooth to your smartphone, it allows full control through the free Vibio app, where you can set multiple alarms, adjust vibration intensity, and receive call or message alerts. Compact, rechargeable, and travel-friendly, Vibio ensures you start every day on time — even if your phone battery runs flat overnight.

### **FEATURES**

- Powerful vibration alarm ideal for deaf, hard-of-hearing, and heavy sleepers
- Wireless Bluetooth connection to iOS and Android devices via the Vibio app
- Create multiple alarms and customise vibration intensity and patterns
- Smart snooze strap – simply pull to delay your alarm without reaching for your phone
- Operates independently even if your phone battery is flat overnight
- Rechargeable lithium battery with USB charging
- Compact, portable design perfect for home or travel use
- NDIS compliant and suitable for accessible alarm needs

### **SPECIFICATIONS**

- **Connectivity:** Bluetooth (compatible with iOS and Android devices)
- **Power Supply:** Rechargeable lithium battery via USB
- **Alarm Type:** Vibration with adjustable intensity
- **Controls:** Through Vibio mobile app

- **Smart Snooze:** Pull-strap activation
- **Operation:** Works independently of phone once alarms are set
- **Use:** Home or travel
- **Compliance:** NDIS approved product

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriiffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)