



BELLMAN VISIT BABY CRY TRANSMITTER



SKU: BE1491

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Visit Baby Cry Transmitter** detects your baby's crying and sends an instant wireless signal to your Bellman Visit receiver, alerting you through vibration, sound, or flashing light. Designed for convenience and peace of mind, it helps ensure you never miss when your baby needs attention, whether during the day or night.

FEATURES

- Wireless baby monitor that alerts you to your baby's crying
- Adjustable sensitivity to suit individual noise and room conditions
- Activates Bellman Visit receivers to provide vibration, light, or sound alerts
- Simple upright placement near the baby for optimal performance
- Ideal for parents who are deaf, hard of hearing, or heavy sleepers
- Can be paired with the Bellman Visit Pager for daytime use or Bed Shaker for nighttime alerts

SPECIFICATIONS

- **Weight:** 0.2 kg
- **Dimensions:** 7.5 × 4.7 × 11 cm
- **Compatibility:** Requires at least one Bellman Visit receiver to operate
- **Power:** Battery operated
- **Placement:** Upright near the baby's cot or crib

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au