




## BELLMAN VISIT DOOR ENTRY (CONTACT) MAT



**SKU:** BE9026

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The **Bellman Visit Door Entry (Contact) Mat** is a practical accessory that helps monitor movement in and out of key areas. When someone steps on the mat, it triggers a connected Bellman Visit transmitter to send an alert to your receivers. Ideal for caregivers, it can be placed beside a bed, doorway, or hallway to provide discreet notifications of activity, enhancing safety and awareness within the home or care environment.

### FEATURES

- Detects movement when stepped on and sends alerts through the Bellman Visit system
- Ideal for use beside a bed or at doorways to monitor entry and exit
- Durable, non-slip surface suitable for daily use in home or care settings
- Provides caregivers with reliable, real-time notifications of movement
- Easy to install and reposition as needed

### SPECIFICATIONS

- **Function:** Movement detection mat for Bellman Visit alerting system
- **Usage:** Triggers alerts when pressure is applied
- **Placement:** Suitable for bedsides, entrances, or exits
- **Compatibility:** Requires Bellman Visit Telephone Transmitter (BE1431)
- **Connection Type:** Wired connection to compatible transmitter

For further information go to **[medimart.com.au](https://www.medimart.com.au)** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)