



### **BELLMAN VISIT DOOR TRANSMITTER**

\$220.00 Original price was: \$220.00.<u>\$198.00</u>Current price is: \$198.00.



**SKU:** BE1411

### PRODUCT DESCRIPTION

# **PURPOSE & BENEFITS**

The **Bellman Visit Door Transmitter** is a smart, wireless alert device that integrates with the Bellman Visit system to notify you when someone rings the doorbell or uses the intercom. Designed for individuals who are deaf or hard of hearing, it detects a wide range of doorbell sounds and sends signals to your Visit receivers, ensuring you never miss a visitor—no matter where you are in your home.

### **FEATURES**

- Automatically detects most doorbell and intercom sounds
- Can be programmed to recognize custom doorbell tones
- Activates Visit receivers with sound, light, or vibration alerts
- Multiple activation methods: internal microphone, built-in electromagnetic detector, external microphone, or direct galvanic connection
- Dual microphone technology for reliable sound detection
- Long-range, robust radio signal covers your entire home
- Battery powered and virtually maintenance-free

# **SPECIFICATIONS**

- Function: Wireless doorbell/intercom detector for Visit alert system
- **Activation Methods:** Internal microphone, electromagnetic detector, external microphone, or galvanic connection
- Power Supply: Standard AA batteries



- Compatibility: Works with all Bellman Visit receivers
- Coverage: Whole-home radio link range
- **Setup:** Simple installation, maintenance-free operation