

BELLMAN VISIT DOOR TRANSMITTER



SKU: BE1411

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Visit Door Transmitter** is a smart, wireless alert device that integrates with the Bellman Visit system to notify you when someone rings the doorbell or uses the intercom. Designed for individuals who are deaf or hard of hearing, it detects a wide range of doorbell sounds and sends signals to your Visit receivers, ensuring you never miss a visitor—no matter where you are in your home.

FEATURES

- Automatically detects most doorbell and intercom sounds
- Can be programmed to recognize custom doorbell tones
- Activates Visit receivers with sound, light, or vibration alerts
- Multiple activation methods: internal microphone, built-in electromagnetic detector, external microphone, or direct galvanic connection
- Dual microphone technology for reliable sound detection
- Long-range, robust radio signal covers your entire home
- Battery powered and virtually maintenance-free

SPECIFICATIONS

- **Function:** Wireless doorbell/intercom detector for Visit alert system
- **Activation Methods:** Internal microphone, electromagnetic detector, external microphone, or galvanic connection
- **Power Supply:** Standard AA batteries

- **Compatibility:** Works with all Bellman Visit receivers
- **Coverage:** Whole-home radio link range
- **Setup:** Simple installation, maintenance-free operation

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au