



### **BELLMAN VISIT PAGER CHARGER**

\$190.00 Original price was: \$190.00.<u>\$171.00</u>Current price is: \$171.00.



**SKU: BE1260** 

## **PRODUCT DESCRIPTION**

#### **PURPOSE & BENEFITS**

The **Bellman Visit Pager Charger** keeps your Visit Pager powered and ready for use while also providing overnight alert functionality. When paired with up to two bed shakers, it ensures you are awakened by vibrations if any Bellman Visit transmitters are activated during the night, offering peace of mind and continuous monitoring while you sleep.

#### **FEATURES**

- Recharges the Bellman Visit Pager battery overnight
- Connects up to two bed shakers for vibration alerts while sleeping
- Automatically triggers bed shaker alerts when Visit transmitters are activated
- Compact, durable, and easy to use for bedside placement
- Ideal for individuals with hearing loss who need reliable night-time alerts

# **SPECIFICATIONS**

• Compatibility: Bellman Visit Pager and Visit transmitters

• Connections: Supports up to two bed shakers

• Function: Pager charging and vibration alert relay

• Power Source: Mains powered

• **Use:** Night-time charging and alerting station