


BELLMAN VISIT PAGER CHARGER



SKU: BE1260



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Visit Pager Charger** keeps your Visit Pager powered and ready for use while also providing overnight alert functionality. When paired with up to two bed shakers, it ensures you are awakened by vibrations if any Bellman Visit transmitters are activated during the night, offering peace of mind and continuous monitoring while you sleep.

FEATURES

- Recharges the Bellman Visit Pager battery overnight
- Connects up to two bed shakers for vibration alerts while sleeping
- Automatically triggers bed shaker alerts when Visit transmitters are activated
- Compact, durable, and easy to use for bedside placement
- Ideal for individuals with hearing loss who need reliable night-time alerts

SPECIFICATIONS

- **Compatibility:** Bellman Visit Pager and Visit transmitters
- **Connections:** Supports up to two bed shakers
- **Function:** Pager charging and vibration alert relay
- **Power Source:** Mains powered
- **Use:** Night-time charging and alerting station

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au