

BELLMAN VISIT PORTABLE RECEIVER



SKU: BE1033

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Visit Portable Receiver** keeps you informed of important alerts throughout your home by providing loud sound signals and bright visual notifications. Ideal for individuals who are deaf, hard of hearing, or need portable alerting, it connects wirelessly to Bellman Visit transmitters to indicate events such as doorbells, phone calls, smoke alarms, or baby monitors. Its lightweight, battery-powered design allows for easy placement anywhere in your home or even outdoors on the balcony or porch.

FEATURES

- Delivers a loud sound signal up to 93 dB for reliable notifications
- Uses coloured lights and distinct sounds to identify each alert source
- Wireless and portable design allows easy movement around the home
- Long-lasting battery life of up to 3 years on standard batteries
- Adjustable radio key supports multiple systems operating in the same area
- Includes wall mount bracket for flexible installation options

SPECIFICATIONS

- **Alert Types:** Audible tones up to 93 dB and coloured light indicators
- **Compatibility:** Works with all Bellman Visit transmitters
- **Power Supply:** Standard replaceable batteries (up to 3 years lifespan)
- **Mounting:** Freestanding or wall-mounted using included bracket
- **Wireless Range:** Typical coverage within standard home environments

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au