

BELLMAN VISIT PUSH BUTTON TRANSMITTER



SKU: BE1240

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Visit Push Button Transmitter** functions as a versatile wireless doorbell or portable personal alarm. When pressed, it sends a signal to your Bellman Visit receivers, alerting you through flashes, sounds, or vibrations. Ideal for use in home or care environments, it ensures you never miss a visitor or call for assistance. Its reliable radio link provides dependable performance across your home, supporting independent living and peace of mind.

FEATURES

- Functions as a wireless doorbell or portable call button
- Sends alerts to Visit receivers via flashes, sound, or vibration
- Customisable signal options for different doors or alert points
- Adjustable radio key allows multiple systems to operate in one area
- Highly reliable radio link ensures consistent transmission
- Battery life of up to 2 years on a standard battery

SPECIFICATIONS

- **Function:** Wireless transmitter for the Bellman Visit alerting system
- **Operation:** Push button activation
- **Transmission Range:** Dependent on environment (may be reduced by brick, concrete, or steel walls)
- **Compatibility:** Works with all Bellman Visit receivers
- **Power Supply:** Standard replaceable battery (up to 2-year lifespan)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au